

# Basic concept

- Dual-view SPIM (like diSPIM) but with horizontal objective geometry
  - Sample is c. elegans worms in vertical tube with thin transparent walls and outer water chamber (customer's design)
  - Main application is GCaMP-type functional imaging
- Use Nikon 40x/0.8 W objectives for light-gathering and axial resolution even though intentionally will undersample laterally for speed
- Two dual-view SPIM (total 4 objectives\*\*) where the pairs are displaced slightly vertically to look at top and bottom half of worm

# Tricky things about the design

- Remote refocus requires alignment adjustment between every pair of lenses (O1-TL1-TL2-O2) for axial spacing, beam tip/tilt, and beam XY
- Remote refocus alignment requirement eliminates possibility of using ASI's standard SPIM head b/c space is at a premium and objectives need to be independently adjusted axially
- Plan: align each remote refocus microscope separately, and subsequently arrange the primary objectives to be co-focused and orthogonal using motorized XY stages and motorized\*\* Z stages
- Using ASI beam scanners for illumination\*\*

# Remote refocus

- Require magnification of 1.333 to remote image
  - O1 = Nikon 40x/0.8 W (5mm EFL)
  - TL1 = ASI C60-Tube-133D dual achromat tube lens (133.3mm EFL)
  - TL2 = ASI C60-Tube-200 achromat tube lens in shorter body (200mm EFL)
  - O2 = Nikon 20x/0.75 air (10mm EFL)
- O3 will be on objective piezo for taking stacks
  - O3 = Olympus 20x/0.8 which is lightweight
  - Plan is for 150um travel piezo = 112um travel in sample space
- Total magnification options quoted per request:
  - With Olympus TL3 =>  $M = 1.333 * 20 \Rightarrow 0.24\mu\text{m}$  pixels,  $\sim 560\mu\text{m}$  square FOV
  - With ASI Tube-100D TL3 =>  $M = 1.333 * 100/9 \Rightarrow 0.44\mu\text{m}$  pixels,  $\sim 1\text{mm}$  square FOV
  - Remember that objective's designed FOV is only  $440\mu\text{m}$  square =  $625\mu\text{m}$  diameter => expect vignetting and increased aberrations outside of that

# Light sheet illumination


- To place MEMS at focal point of 200mm lens:
  - Scan-SL without scan lens
  - modified C60-SPACER-ADJ for axial positioning and adapted to Scan-SL
- max beam NA
  - $1.2\text{mm} \cdot \cos(22.5)/2$  MEMS size
  - $/200 \cdot 133/5 = 0.073$  max NA
  - $0.64 \cdot 0.488\mu\text{m} \cdot 1.333 / .073^2 = 78\mu\text{m}$  confocal length => plenty of NA
    - NA 0.029 gives 495 $\mu\text{m}$  confocal length
    - => consider 7.5mm collimator
- want ~500 $\mu\text{m}$  scan range to cover objective FOV
- Travel range at sample is OK
  - +/- 0.14 radian MEMS range (optical)
  - 200mm TL2
  - 133mm TL1
  - 5mm O1
  - => 1.05mm FOV or about 2x required

# Imaging speed

- Application calls for rapid volumetric imaging
- Camera readout is fundamental speed limit with readout requiring  $\sim 5\mu\text{m}/\text{row}$
- Large pixels mostly to have smaller images for faster readout
- Hamamatsu Fusion selected for good speed and excellent all-around performance and versatility
- 1.33x mag to remote image, additional mag of 20x or 11x to sensor => total mag will be 26.7x or 14.8x
- Objective-limited FOV of  $\sim 440\mu\text{m}$  square =  $\sim 625\mu\text{m}$  diagonal
- $440\mu\text{m}$  square =  $\sim 1800$  pixels square at 26.7x;  $\sim 1000$  pixels square at 14.8x

# Layout

**Side view**




O1/TL1 spacing  
 133D will have ~111mm free focal length towards collimated space  
 13 inside objective  
 8 bushing  
 0-4 bushing gap (adjustable)  
 14? XY adjuster  
 60 Cube-III w/ mirror  
 ~14 extra space => C60-EXT-15


TL1/TL2 spacing  
 133D will have 60mm free focal length towards focus space  
 200 in 100 body will have ~163mm free focal length towards focus space  
 37 AFS body minimum  
 0-8 AFS travel (adjustable)  
 60 Cube-III with mirror  
 60 Cube-III with mirror  
 ~62 extra space => EXT-37 + EXT-25

TL2/O2 spacing  
 200 in 100 body will have ~190mm free focal length towards collimated space  
 10 inside objective  
 8 bushing  
 0-4 bushing gap (adjustable)  
 14? XY adjuster  
 60 Cube-III w/ mirror  
 60 Cube-III w/ dichroic  
 ~36 extra space => EXT-25 + EXT-10

**Top view**



TL2/LS spacing  
 200 in 100 body will have ~190mm free focal length towards collimated space  
 60 Cube-III  
 87 Spacer-Adj body  
 0-30 travel  
 30 minimum Scan-SL body  
 13 extra space



# Estimated mass

- Estimate 10 kg for assembly to be raised/lowered
- LS-50 limit is 5 kg each
- => right near limit but should be OK if we take steps to e.g. trim breadboard

<u>Item</u>	<u>Mass [g]</u>	<u>Qty</u>	<u>Total</u>	<u>Comments</u>
breadboard	3000	1	3000	10 g/mm with 300mm width, could trim
Cube-III	400	5	2000	including mirror, includes spacer-ADJ cube
spacer-ADJ	600	1	600	
scanner	700	1	700	
AFS	600	1	600	
O1	120	1	120	
O2	250	1	250	
OBLPA	150	2	300	stand-in for XY positioner
Ext-75	225	2	450	Stand-in for tubes
cube jack	75	2	150	
Tube 133D	420	2	840	
mounting	300	3	900	guess
total			9910	



# Remote refocus alignment procedure Part 1

- Requires collimated light source (beam size options ~5mm and ~20mm) and shear plate
  - not included in quote but can be provided by ASI if desired
  - Procedure can be tedious but it is algorithmic
- O1/TL1
  - Shine collimated light through TL1 towards O1 and check back-reflection
  - Use O1 bushing to get collimation
  - Use M1 and O1's XY to get non-deflection when O1 added/removed
- TL1/TL2
  - Shine collimated light through TL1 towards TL2 and check back-reflection
  - Use AFS-800 to get collimation
  - Use M2 and M3 to get non-deflection when TL2 added/removed


# Remote refocus alignment procedure Part 2

- TL2/O2
  - Shine collimated light through TL2 towards O2 and check back-reflection
  - Use O2 bushing to get collimation
  - Use M4 and O2's XY to get non-deflection when O2 added/removed
  - (do this with dichroic in place, but dichroic doesn't need to be aligned yet)
- O2/O3
  - Shine collimated light through TL1 through system and check back-reflection
  - Move O3 in XYZ to find focus of that light


# Objective alignment procedure

- Adjust height of objectives and lateral positions by eye, then turn on cameras with dye solution in chamber
- Use XY and Z stages to adjust 3D position of each objective

# Remote refocus aberrations seem tolerable



# Spot diagram of remote refocus



Aberrations are quite minimal, just start to vignette at objective full field of 625um diameter

# New parts for ASI

- APZOBJ-150-S
  - Extension of existing product
  - Parts are already being made and should be ready to go shortly
- Tube-133D
  - Optic design done, will be ~1 month from order to make parts
- XY adjuster
  - Require adjuster akin to Thorlabs CXY1 but compatible with ASI's bushings
  - Will probably be designed for other purposes anyway
  - Have sketch, need to do full mechanical design and machine
  - Ready 6-8 weeks from when we get serious about it (an order makes it serious)